Medical Imaging

Medical Imaging

- IBM Visual Insights is a researcher tool to do rapid labeling and deep learning classification of medical images.
- VI runs on Intel x86 systems and IBM POWER8/9 Systems
- It is being used actively across multiple disciplines. We have several MVP use cases developed.
- See Medical Imaging Guide.


The IBM medical imaging solutions can provide speedup of image analysis to meet clinical, medical, or academic researchers' needs. (Fig 1).

Fig 1. Examples of medical imaging applications

	Dermatology	Diagnostic Radiology	Digital Pathology
Data Characteristics	Photos: RGB, 2D	X-rays, CT scan, ultrasound, MRI	High res digital slides, 2D, 3D
	Malahoma Atgicel Berign		0.00 mg
Imaging Challenge	Identifying malignant mole with mobile phone images	Detecting lung nodules on 3D CT scans	Detecting and classifying blood cell subtypes in blood samples
Desired Outcomes	Prediction of malignancy in shorter time frame	Detection of lung cancer with higher accuracy	Diagnosis of blood- based diseases efficiently
Techniques also applicable to	Ophthalmology Dentistry	Breast imaging Cardiothoracic Cardiovascular Chest Endovascular Gastrointestinal Genitourinary Head and neck Musculoskeletal Neuroradiology Nuclear Radiation oncology Dentistry	Chemical Clinical Cytopathology Dermatopathology Forensic Genetic Hematology Immunopathology Medical microbiology Molecular Neuropathology


Visual Insights: "Point-and-Click" Al for Images & Video

Label Image or Video Data






Patient without DR Patient with DR



Auto-Train AI Model

Training / NWPU-Breast Tumor (magnification 40×)_model




Package & Deploy AI Model










Breast Cancer Classification





nd models are created from a nodel can be validated, expor	prepared this sets. ried, and deployed for production to	for Graphics Processing Units (GPU).		el		
taper	f model					
ge classification set IMPU dress Jump I set Debut theig elect	inspringtin 40%	Created 12/16/2019, 1:13:07 FM By admin				
88% Accuracy	Model hyperparan Max iteration @ 1500 Test interval @ 20 Weight decay @ 0.0005	Test iteration @ 100 Learning rate @ 0,001 Ratio @ 0,8				
86% Progration)	S8% Recuil	Loss VS Iterati	on		
Programme		88%	-		Train Loss Train Loss Train Loss	
ategory)	88% Recuil	-		Train Loss Train Loss Train Loss	
Alegory -ductal_carchiema St-inuges -filtrosdenema	Precision	SB% Recall	-	Sandon J. J. J. J.	Train Loss Train Loss Train Loss	
Inggry Inductal carchinema Shimupin Inductal carchinema Shimupin Inductal carchinema Li magin Inductal selections	Precipion 0.895	88% Recall 0.837	-	Sandon J. J. J. J.	Train Loss Train Loss Train Loss	
86%, Precision of the Company of the	Precision 0.895 0.109	8856, Recall 8.837 0.909	-	Sandon J. J. J. J.	Train Loss Train Loss Train Loss	







Diabetic Retinopathy Detection


Use Case Description: Diabetic retinopathy (DR) is the leading cause of blindness. DR will grow from 126 million in 2010 to 191 million by 2030. Comprehensive and automated method of DR screening has long been recognized.


Solution: Create an automatic image classification solution to assist medical diagnosis and research etc.

Dataset: https://www.kaggle.com/c/diabetic-retinopathy-detection

35,000 + images of various classes

- 0 No DR
- 1 Mild
- 2 Moderate
- 3 Severe
- 4 Proliferative DR





Diabetic Retinopathy Detection

Step 1
Prepare Data

Step 2
Generate Dataset

Step 3

Data Augmentation

Step 4
Train Models

Step 5

Deploy Models

Step 6
Inferencing

Catalog large number of data files through metadata management. It is almost impossible to manage large number of files manually.

Automatically **select any dataset** and subset for multiple training iterations.


Create dataset with **balanced** categories to improve model accuracy. Simple UI to create new image data.

Automatic model training using **pre-trained model** and **GPU** to scale compute.

One click model API end-point deployment for any application.

Infuse AI with online (mobile/web app), batch, and real-time streaming process.
























Scale up to Fight Covid-19

Covid-19 Research and Diagnosis



Training and Inferencing

offset sex age finding survival intubated intubation present went icu needed supplemental O2 extubated temperature pO2 saturation leukocyte count neutrophil count lymphocyte count view modality date location folder filename doi url license clinical notes

other notes


Covid-19 Research and Diagnosis

Answer Key Questions from Scientific Literature



- What is known about transmission, incubation, and environmental stability?
- What do we know about COVID-19 risk factors?
- What do we know about virus genetics, origin, and evolution?
- What has been published about medical care?

Read COVID-19 X-ray or CT image



While PCR tests offer many advantages they are physical things that require shipping the test or the sample. X-ray machines can be plugged in to screen patients as long as they have electricity.

Al tools can help general practitioners to triage and treat patients. Companies are developing Al tools and deploying them at hospitals <u>Wired 2020</u>.

Source: <u>Kaggle</u> Source: <u>IEEE</u>

Open Issues

- Al Golden Age: Big Data, Big Compute, Deep Learning
- State-of-the-art deep neural network architectures
 - Supervised pre-training/domain-specific fine-tuning
 - Unsupervised learning GAN (Generative adversarial network)
 - Data + knowledge
- Open medical datasets: From ImageNet to ??
 - Data Labeling and Annotation
 - Data Privacy Federated Learning
- High Performance Computing in FPGA, GPU, or CPU (x86 or P9)
 - Benchmarking MLPerf Medical Imaging
 - Model size reduction
 - Deployment on Edge Devices

Questions?

Twitter: @a9zhang

Email the speaker: <u>andrew.zhang@ibm.com</u>

Call to Action

- Contact Ganesan or AICOC (aicoc@us.ibm.com) for guided lab exercises
- Collaborate to share data and models for both DR and X-Ray lung Covid-19
- Looking for PhD students with medical background to develop advanced research on IBM super-computers
- Post slides on my twitter and Linkedin account
- Twitter: @a9zhang
- Email the speaker: andrew.zhang@ibm.com



— HEALTHCARE Issued on: March 23, 2020

STATEMENTS & RELEASES

Call to Action to the Tech Community on New Machine Readable COVID-19 Dataset

— HEALTHCARE Issued on: March 16, 2020





References

- 1. Rapid AI Development Cycle for the Coronavirus (COVID-19) Pandemic: Initial Results for Automated Detection & Patient Monitoring using Deep Learning CT Image Analysis
- 2. COVID-19 Screening on Chest X-ray Images Using Deep Learning based Anomaly Detection
- 3. <u>Vendors debut AI X-ray system for COVID-19, give it away for free</u>
- 4. How Does COVID-19 Appear in the Lungs?
- 5. <u>COVID-19</u>

